


## eLAB Control Unit

Weight and dimensions		
	Single	~ 35 kg
	Twin	~ 50 kg
	Single / Twin	437 x 850 x 500 mm
Hardware		
	Housing	Stainless Steel, AISI 304 (SF4 according to ASME BPE)
	Display	Capacitive Touch Screen, 15", glass
	Status indicator light	Green: OK Orange: warning Red: alarm
Software		
	eSCADA Basic	Included
	eSCADA Advanced	Optional
	CFR21.11	Optional
	Remote access	Optional
Documentation package		
	GMP	Optional
Built-in peristaltic pumps (4x)		
	3x fixed speed: 5 – 200 rpm	1x variable speed: 5 – 600 rpm
	* External pumps optional (4x)	
	Addition bottles	Optional
	Tubing	Optional
Process Control Sensors		
	Temperature	Pt100: 10 – 60 °C
	Dissolved Oxygen	Optical sensor: 0 – 100 % (ODO Cap membrane as spare part)
	pH Sensor	Glass, pre-pressurized reference system: 2 – 12 pH
	Level/Foam	Hygienic. Sliding connection
	Redox	Optional
	Conductivity	Optional
	Total Cell Density (TCD)	Optional
	Viable Cell Density (VCD)	Optional
	Biochemical parameters	Optional
	pCO <sub>2</sub>	Optional
	Weight control	Optional
	Digital/Analogic auxiliary input	Up to 4 4 – 20 mA / RS485 Modbus protocol
Motor		
	Agitator	Magnetic agitation
	Rotation Speed Motor	1 L: 0 – 2000 rpm; 2 L: 0 – 1800 rpm; 5 L: 0 – 1500 rpm;
	Impellers	Microbial: 2x or 3x Rushton Cell culture: 1x or 2x Marine
	*Customized impellers available	
Temperature Module		
	Single wall	Electrical heating (blanket) and cooling loop with external chiller
	Double wall	Thermostat system with external chiller for cooling and electrical resistance for heating
Aeration Module		
	Microbial	Three-gas mixing: Air, O <sub>2</sub> , (N <sub>2</sub> optional) Three flow meters through sparger by MFCs Gas mix as aeration strategy: AIR and O <sub>2</sub> Up to 2 vvm volume flow for each gas
	Cell Culture	Four-gas mixing: Air, O <sub>2</sub> , CO <sub>2</sub> , N <sub>2</sub> by MFCs Four flow meters through sparger and overlay Up to 1 vvm volume flow for each gas
Utilities		
	Power supply	230 V (EU)
	Gases	1 bar
	Water	Water supply pressure from external chiller, 1 bar



TECNIC's bioreactors are available with multiple configuration possibilities for microbial and cell culture applications. All the equipment is designed and manufactured in our facilities in Riudarenes (Girona), Spain.

Our bioreactors are available for laboratory, pilot and production applications, fitting in all processes from small scale to final product production.



### eSCADA

Our software, built on the AVEVA (formerly Wonderware) System Platform architecture, is designed to facilitate the recording and control of all process parameters in an intuitive and easy-to-use way.

eSCADA is integrated into all of our equipment, which means that there is full interconnectivity between them. No external software is required to control any part of the process. The system can also be remotely controlled from an external PC or tablet.

**eSCADA Basic** is the standard software that runs the equipment, where the process can be controlled, visualized and graphed, as well as export results in an efficient way.

**eSCADA Advanced** (optional) is the upgraded version for process automation that allows recipes management and control loops. Both versions can include GMP and CFR 21 Part 11 modules.

### pH

pH regulation is performed by a loop with acid and base pumps (CO<sub>2</sub> is also available for cell culture) based on a setpoint. Digital technology is used.

### pO<sub>2</sub>

pO<sub>2</sub> regulation is performed by cascade. Up to 4 actuators: air, oxygen, stirring and substrate. The user can select its order depending on the culture. Digital technology is used.

### °C

Temperature is regulated by a thermostat system with external chiller for cooling and electrical resistance for heating (heating blanket and cooling finger for single wall vessels are available).

BIOREACTORS CONFIGURATION	eLAB		ePILOT		ePROD	
	RB-M-0001/2/5	RB-C-0001/2/5	RB-M-0010/50	RB-C-0010/50	RB-M-0100/2000	RB-C-0100/2000
<b>1.VESSEL</b>	<ul style="list-style-type: none"> <li>Single/Twin Glass: <input checked="" type="radio"/> STD <input type="radio"/> OPT</li> <li>Stainless Steel: <input type="radio"/></li> <li>Single Use: <input type="radio"/></li> <li>Working Volume (L)*: 1-2-5 (eLAB), 10-30-50 (ePILOT), 100-200-500-1000-2000 (ePROD)</li> <li>Weight Control: <input type="radio"/></li> <li>Single/Double Wall: <input checked="" type="radio"/> STD <input type="radio"/> OPT</li> <li>Lighting: <input type="radio"/></li> </ul>					
<b>2.AGITATION</b>	<ul style="list-style-type: none"> <li>Rushton Impeller: <input checked="" type="radio"/></li> <li>Marine Impeller: <input type="radio"/></li> <li>Customized: <input type="radio"/></li> <li>Deflectors (4): <input type="radio"/></li> <li>Magnetic Agitator: <input checked="" type="radio"/></li> </ul>					
<b>3.ADDITION</b>	<ul style="list-style-type: none"> <li>Fixed Speed (3x) (ml/min): <input checked="" type="radio"/></li> <li>Variable Speed (1x) (ml/min): <input checked="" type="radio"/> Substrate / Media</li> </ul>					
<b>4.MEASUREMENT</b>	<ul style="list-style-type: none"> <li>pH: <input checked="" type="radio"/></li> <li>pO<sub>2</sub>: <input checked="" type="radio"/></li> <li>Temperature: <input checked="" type="radio"/></li> <li>Level: <input type="radio"/></li> <li>Foam: <input checked="" type="radio"/></li> <li>Conductivity: <input type="radio"/></li> <li>ORP [-1500 - 1500]mV: <input type="radio"/></li> <li>Cell Density (OD, viability): <input type="radio"/></li> <li>Biochemical Parameters (glucose...): <input type="radio"/></li> <li>pCO<sub>2</sub> (10 - 1000) mbar: <input type="radio"/></li> </ul>					
<b>5.GASSING</b>	<ul style="list-style-type: none"> <li>Air Flow (lpm): <input checked="" type="radio"/></li> <li>O<sub>2</sub> Flow (lpm): <input checked="" type="radio"/></li> <li>N<sub>2</sub> Flow (lpm): <input type="radio"/></li> <li>CO<sub>2</sub> Flow (lpm): <input checked="" type="radio"/></li> <li>Gas Mix (%): <input checked="" type="radio"/></li> <li>Gas Inlet: <input checked="" type="radio"/> Sparger <input type="radio"/> Sparger / Overlay</li> </ul>					
<b>6.SOFTWARE</b>	<ul style="list-style-type: none"> <li>eSCADA Basic: <input checked="" type="radio"/></li> <li>eSCADA Advanced: <input type="radio"/></li> <li>CFR21.11: <input type="radio"/></li> </ul>					
<b>7.DOCUMENTATION</b>	<ul style="list-style-type: none"> <li>GMP: <input type="radio"/></li> </ul>					



15" Capacitive touch screen

Standard version: 4 peristaltic pumps

Twin version: 8 peristaltic pumps (Watson Marlow, 114DV; flows up to 510 mL/min; no maintenance; reversible direction)

pH, pO<sub>2</sub> and foam/level sensors by default

Light indicator for a quick status view

Magnetic agitation

Gas supply systems for sparger and overlay controlled by four flow meters: air, oxygen, nitrogen and carbon dioxide. Gases distinguish between cell culture, aerobic and anaerobic microbial

For anaerobic processes, air and oxygen must be changed for nitrogen



### eLAB Vessel

Our benchtop bioreactors are made of borosilicate glass with working volumes of 1, 2 and 5 L. They are available in single or double wall

Optional Stainless Steel vessel for 2 and 5 L which may include Cleaning in Place (CIP) cycle module

Glass Material	Standard Vessel: Borosilicate Glass Wetted parts: Stainless Steel AISI 316L (SF1 according to ASME BPE) Other parts: Stainless Steel AISI 304L (SF4 according to ASME BPE)
Size (working volume)	1 L (0.35 - 1 L) 2 L (0.6 - 2 L) 5 L (0.8 - 5 L)
Dimensions (diameter x height) <i>*With structure; without condenser</i>	1 L - 210 x 300 mm 2 L - 235 x 400 mm 5 L - 260 x 500 mm
Sterilization	Autoclave
Baffles	4x integrated baffles. Optional
Sparger	Ring sparger <i>*Microsparger optional</i>
Stainless Steel 316L	Optional (2 or 5 L)
Cleaning-in-place (CIP)	Optional