


TECNIC

BIOPROCESS EQUIPMENT MANUFACTURING

eLAB
ePILOT

TANGENTIAL FLOW FILTRATION



TECNIC's fully automatic tangential flow filtration equipment has a control unit for processing and a vessel (optional) for media concentration or diafiltration. This vessel is available in stainless steel for reusable or in polycarbonate for single-use applications. The full design and production process is performed at our facilities in Riudarenes (Girona), Spain.

Available for laboratory, pilot and production applications, our equipment fit in all processes from small scale to final product production.


eSCADA

Based on a Supervising Control And Data Acquisition (SCADA) type system and built of the Wonderware System Platform architecture, this software allows the management of basic process control parameters up to advanced recipes and report management.

By means of loops, this software allows to control variables such as pH (acid, base, CO₂), pO₂ (up to 4 cascades including gases and agitation), temperature and antifoaming, as well as the prediction of cell behaviour at metabolic level with in-situ cell viability sensors

eSCADA BASIC (optional)

The standard software that runs the equipment, where the process can be controlled, visualized and graphed, as well as export results in an efficient way.



eSCADA ADVANCED

The upgraded version for process automation that allows recipes management and control loops. Both versions can include GMP and CFR 21 Part 11 modules.



Flow Sensor

eLAB Vessel (optional)
The vessel is made of stainless steel and is available in 5 and 10 liters. Its bottom has been designed to optimize the minimum volume when harvesting. As the temperature is an important process parameter, double wall tanks are also available. Double wall temperature is regulated by a thermostat system with external chiller for cooling.

Stainless steel vessels are self-cleaning as 360° spray balls are integrated so Cleaning in Place (CIP) module is included in the TFF system.

Polycarbonate vessels are also available if the process requires Single-Use containers.



Automatic Control Valves

Pressure Sensor

Holder and cassettes

10" Capacitive touch screen with the software eSCADA

Two Peristaltic Pumps: one has fixed speed and another one variable

Light indicator for a quick status view (green, orange and red)

Transmembrane Pressure (TMP) Control

This process parameter is the average applied pressure from the feed on the permeate side of the membrane.

$$TMP [bar] = \left(\frac{P_{feed} + P_{retentate}}{2} \right) - P_{permeate}$$

This eLAB crossflow system is fully automated, so TMP can be defined previously and will be controlled by acting on the retentate valve throughout the process.

Permeate Flux Constant

Filtrate flux is the filtrate flow rate normalized for the area of the porous membrane through which it passes. By controlling permeate valve, the system will automatically regulate all the process.

ΔP Control

ΔP is a process parameter that is defined as the difference between feed inlet and retentate outlet pressure of the membrane.

$$\Delta P [bar] = P_{feed} - P_{retentate}$$

This value is controlled automatically once the set point is set by the recirculation pump.

eLAB Control Unit

Hardware	Housing Display	Stainless Steel, AISI 304L (SF4 according to ASME BPE) Capacitive touch Screen, 10", glass
Software	eSCADA Basic eSCADA Advanced CFR 21.11 Remote access	Included Optional Optional Optional
Documentation Package	GMP	Optional
Built-in peristaltic pumps	Fixed speed <i>*Feed, permeate, diafiltration</i> Variable speed <i>*Feed, permeate, diafiltration</i> Recirculation pump <i>*External pumps optional</i>	1x fixed speed: 90 rpm 1x variable speed: 100 rpm 1x sanitary four-piston pump: up to 1400 L/h
Process Control Sensors	Level (vessel) Flow sensor Pressure Temperature (vessel) Load Cells (vessel) pH Sensor (feed) Conductivity (feed) Magnetic agitation (vessel)	Guided radar sensor 2x flow sensor for feed/permeate 3x pressure sensor for feed/retentate/permeate Optional (Pt100 / Pt1000) Optional Optional Optional Optional
Temperature module	Single wall Double wall	No temperature control. Visualization Thermostat system with external chiller for cooling
Tangential Flow Filtration System	Holder Porous membrane Cut off Filter area Cassettes <i>*Typology depends on product</i> Minimum recirculation volume Maximum inlet pressure (SS piping)	eHOLDER by TECNIC Sartorius & Pall & Millipore compatible Microfiltration / Ultrafiltration / Nanofiltration 0.2 0.45 µm / 2 5 10 30 100 300 kDa 0.1 – 2 m ² 1 to 5 ~ 150 mL 4 bar
Other filtration modules	Ceramic membranes Hollow Fiber	Optional Optional
Utilities	Power supply Water for cooling	230 V Water supply pressure 1 bar

eLAB Vessel (optional)

Material	Stainless Steel	Optional Wetted parts: SS AISI 316L (SF1 according to ASME BPE) Other parts: SS AISI 304L (SF4 according ASME BPE)
Volume	5 L <i>*Customized volume available</i> 100 L	
Bottom	Conical, 45°	
Single wall	Without temperature control loop	
Double wall	Optional	
Lightning	Optional	
Sight glass	Optional	
Single Use	Optional	

TANGENTIAL FLOW FILTRATION CONFIGURATION

*Other customized options available

	MANUAL	AUTO	AUTO
1.VESSEL			
Without	00	00	
Standard 5 / 10L (h:d, 3:1)	01, 05, 10	05, 10	
Customized (> 10L)			Permeate Tank optional (See Buffers datasheet) Agitation (Optional)
Single Wall			
Double Wall			
2.FILTER			
Cassette	MF, UF, NF	#01, 02, 03, 04, 05	
Ceramic			
Hollow Fiber			
3.TEMPERATURE			
Pt100 Sensor			
Pt1000 Sensor			
4.CONTROL			
Level Sensor			
Load Cells			
Flow Sensor (2x)			
Pressure (3x)			
pH Sensor			
Conductivity			
Agitation			
5.SOFTWARE			
eSCADA Basic			
eSCADA Advanced			
CFR21.11			
Remote Access			
6.DOCUMENTATION			
GMP			